

Sales Event for NEW products



Get a special discount off
To 31st Mar. 2021. (Based on ordered date.)

NEW Product #1.

Fire storage cabinet (SC3-F)

SC3-xxF cabinet has all the requirements that a fire safety storage cabinet needs according to the principles of flammable liquids and material storage

 [YouTube](#) “Real fire simulation test”



Fire resistance performance of the Lab Companion's fire safety cabinet was verified in real fire simulation.



- Three kinds of fire resistance performance type - 30, 60, and 90min. (EN 14770-1 certified)
- Various shapes and sizes - Stand or under counter type, available to choose shelves or drawer.
- Optional accessories - Extraction module and its dedicated filters.

NEW Product #2.

Electrophoresis System & Power Supply (EP & EPS)

We suggest a combination use of EP&EPS designed to solve some inconvenient concerns of Electrophoresis experiment, it bring a quicker and easier your experiment.



Observe the gel conveniently without separating it from tray



Multichannel pipettes can be used (31, 16 well combs)




Easy to measure the thickness of gel due to engraved level marker



Equipped with a spirit level Anti-slip rubber feet (EP-10, EP-18 models)



Anti-slip rubber fee



Compatible with other companies' products and can be connected to up to 4 devices at the same time

- Large quantity of samples can be loaded quickly – Available to use Multichannel pipettes for 31 and 16 well combs.
- No need to separate a gel from the tray – Provide UV transmitting gel tray.
- Minimized experiment variation – Level marker and Spirit level for making a constant and flat gel.
- Safety features to reduce mistakes – Clear observation of the direction of progress, and equipped with anti-slip rubber feet.
- Possible to connect an Electrophoresis system up to 4 devices to EPS at the same time.

NEW Product #3.

Safety Waste Disposal System (SF-14C/P, WT-20C/P)

Combination safety funnel and waste container prevent unexpected accidents that may occur by indiscriminate disposal of used chemicals.

 [YouTube](#) “Safety funnel operation video”



Level indicator operation of Lab Companion Safety Funnel for waste disposal.

- General and Anti-static type.
- Level indicator of funnel - Easily check the waste liquid level in the container.
- Various colored lids and classified color labels of funnel.
- Prevent mixing of waste liquids that have different chemical properties.
- Optional accessory - Carbon filter to remove malodor and harmful gas.