

 **Lab Companion**

Quick & Easy to Replace

LC Magnet System

Magnetic Arm Hood

1


Lab Companion New Technology of Arm Hood

What is LC Magnet System?

The new system of the arm hood to replace easily a dome without any assembly thanks to the magnetic on the connection part of the arm and dome.



Lab Companion



2

LC Magnet System

CONVENIENT! QUICK!



- ✓ Magnet type to **replace a dome easily** in 5 seconds
- ✓ No need the additional assembly, **convenient to replace**
- ✓ Sturdy attach with **strong magnet** securing the fume leakage

3

LC Magnet System Makes you more efficient!

Dome type of hood should be selected depending on the properties of gas and the diffuse range of fume. LC Magnet will be a solution to increase space efficiency in case the space is too small to place several arm hoods because the one hood is enough to use for easy dome replacement. Also, it makes to increases the experiment efficiency, no need to use a separate arm hood for each experiment.



Space-saving



Multipurpose



Experimental
Efficiency



Lab Companion